

**Crude
Petrol**

Crude Petrol

Standard ISO 9030:
with 8" conical tube - TU 030
and tube stoppers - PV 156

**Lubricating
oil**

Lubricating oil

Standards:
ASTM D91 with 8" conical tube -TU 030
ASTM D893 with 8" conical tube - TU 030
ASTM D2273 with 8" trace tube - TU 056
ASTM 2711 with 8" conical tube -TU 030

**Crude
Oil**

Crude Oil

Standards:
ASTM D96 with:
8" conical tube -TU 030
or 6" conical tube -TU 050
ASTM D4007 with:
conical tube 8" - TU 030
and tube stoppers - PV 156

Biodiesel

Biodiesel

Standard ASTM D2709:
with 8" trace tube - TU 056

Fuel oil

Fuel oil

Standards:
ASTM D1796 with 8" conical tube - TU 030
ASTM D2709 with 8" trace tube - TU 056
ISO 3734 with 8" conical tube - TU 030

Polishes

Polishes

Standard ASTM D1290:
with 8" conical tube TU 030



Digtor 22 C Series Applications