

PROCEDURE FOR RECYCLING EQUIPMENT (OTHER THAN SPAIN)

Within Spanish territory, Ortoalresa provides an equipment end-of-life recycling system through an integrated waste management system (Directive 2012/19/EU), as a member of the ECOASIMELEC foundation.

However, these waste management systems are only mandatory within Spain. In its commitment to the environment, Ortoalresa proposes a series of recommendations to generate minimal environmental impact and facilitate parts recycling in these territories out of Spain.

- Decontaminate the equipment prior to performing any action and issue the corresponding certificate (this certificate is available on our website, www.ortoalresa.com/descargas).

- In countries with a comprehensive waste management system, contact the managing entities for details of the procedures to follow.

- In countries with no comprehensive waste management system, proceed to destroy/recycle the equipment in accord with the laboratory equipment recycling rules/criteria established in that country.

Any questions or queries, please contact us via e-mail at sat@ortoalresa.com or tel. +34 91 884 40 16